

VisionArray HPV Chip 1.0

REF VA-0001-10 ∇_{Σ} 10 tests

REF VA-0001-50 ∇_{Σ} 50 tests

For the specific detection of 41 Human Papilloma Virus (HPV)-types that have been produced with the help of the VisionArray HPV Primer Kit.

CE

IVD

In vitro diagnostic medical device

according to EU directive 98/79/EC



Trademarks: VisionArray® is a trademark of 42 life sciences GmbH & Co. KG

Manufactured by:

42 life sciences GmbH & Co. KG
Fischkai 1
27572 Bremerhaven/ Germany
Phone: +49 471 4832-500
Fax: +49 471 4832-308
www.42ls.com
Email: info@42ls.com

Distributed by:

ZytoVision GmbH
Fischkai 1
27572 Bremerhaven/ Germany
Phone: +49 471 4832-300
Fax: +49 471 4832-509
www.zytovision.com
Email: info@zytovision.com

Rev.: 1.1

As of: April 01. 2016

Contents

- 1 Product Description 1
 - 1.1 Scope of Application 1
 - 1.2 Product..... 1
- 2 Basic Principles 2
- 3 Safety Precautions and Disposal..... 2
- 4 The VisionArray HPV Chip 1.0..... 3
 - 4.1 Components..... 3
 - 4.2 Storage and Shelf Life 3
 - 4.3 Test Material 3
 - 4.4 Additional Material 3
 - 4.5 Description of the Chip 4
- 5 Interpretation of Results 5
 - 5.1 General..... 5
 - 5.2 Software-Based Evaluation 5
 - 5.3 Specificity of the Tests 6
 - 5.4 Sensitivity of the test 7
 - 5.5 Cutoff 9
 - 5.6 Limits of the Method 9
- 6 Literature 9
- 7 Problems and Possible Causes..... 10
- 8 Explanation of the Symbols 11

1 Product Description

1.1 Scope of Application

The VisionArray HPV Chip 1.0 is intended to be used with the VisionArray Analysis Package for the qualitative detection and genotyping of PCR amplicates of 41 clinically relevant Human Papilloma Virus (HPV) genotypes that have been produced with the help of the VisionArray HPV Primer Kit and the VisionArray Detection Kit.

Infections with HPV are common and a major risk factor for the development of e.g. cervical carcinoma. At present, there are more than 150 different HPV-types described. Depending on their risk to induce cancer, they are divided into Low Risk (LR) and High Risk (HR) types.

The VisionArray HPV Chip 1.0 is designed to detect the following 41 genotypes:

Risk classification of the 41 HPV genotypes on the VisionArray HPV Chip 1.0

High Risk	Probably High Risk	Low Risk
16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59	26, 34, 53, 66, 67, 68a, 68b, 69, 70, 73, 82 (IS39), 82 (MM4)	6, 11, 40, 42, 43, 44, 54, 55, 57, 61, 62, 72, 81 (CP8304), 83 (MM7), 84 (MM8), 90, 91

The HPV types were classified according to the current scientific literature (refer section 6 Literature).

This product is designed for *in vitro* diagnostic use (according to EU directive 98/79/EC). Interpretation of results must be made within the context of the patient's clinical history with respect to further clinical and pathologic data of the patient by a qualified pathologist.

1.2 Product

VA-0001-10: 10 tests

VA-0001-50: 50 tests

2 Basic Principles

DNA fragments with a specific sequence are detected from a pool of DNA fragments on a glass chip with the help of immobilized DNA capture sequences by DNA/DNA-hybridization. For this detection system DNA samples from formalin-fixed, paraffin-embedded tissue or cell samples can be used as raw material. As a first step, the target sequences in these samples have to be amplified and biotinylated in a PCR. The hybridization between the amplified sequences and the complementary DNA captures is performed subsequently. After the hybridization the unspecifically bound DNA is washed away by short stringent wash steps. The specific bound biotinylated sequences are secondary labeled with a Streptavidin-Peroxidase-Conjugate afterwards and visualized by a tetramethylbenzidine (TMB) staining.

3 Safety Precautions and Disposal

- ✓ Read the operation instructions prior to use!
- ✓ Do not use the chips after the expiry date has been reached!
- ✓ Chips should be used in a dust-free setting. Avoid the contamination of the chip surface with dust or other particles!
- ✓ Avoid direct contact with the array field on the chip-surface!
- ✓ Only the labelled side of the slide can be used for hybridization.

4 The VisionArray HPV Chip 1.0

4.1 Components

The following components are included:

Code	Components	Quantity	
		10	50
VA-0001	VisionArray HPV Chip 1.0	10	50
	Instruction manual	1	1

4.2 Storage and Shelf Life

The chips have to be stored in the intact original packing at room temperature. If these storage conditions are followed, the chips are stable, without loss of performance, at least until the expiry date printed on the label.

After opening the original packaging use the chips within one month.

4.3 Test Material

Starting material for this detection system are PCR amplification products that have been produced with the VisionArray HPV Primer Kit (Article-No. VP-0001-50).

The hybridization and detection of the chips has to be performed with the VisionArray Detection Kit (Article-No. VK-0003-50) according to the specifications of the manufacturer.

4.4 Additional Material

The chips have to be analyzed on the VisionArray Analysis Package (Article-No. E-4060-1).

The VisionArray Analysis Package has to contain the VisionArray HPV Chip File 1.0 for a successful scan.

4.5 Description of the Chip

Positioning of the capture sequences on the chip:

1		2	3	4	5	6	7		1
		8	9	10	11	12	13		
14	15	16	17	18	19	20	21	22	23
24	25	26	27	28	29	30	31	32	33
34	35	36	37	38	39	40	41	42	43
11	10	9	8	7	6	5	4	3	2
21	20	19	18	17	16	15	14	13	12
31	30	29	28	27	26	25	24	23	22
		37	36	35	34	33	32		
1		43	42	41	40	39	38		

Allocation position and capture sequence:

No.	Capture sequence	No.	Capture sequence	No.	Capture sequence	No.	Capture sequence
1	Guide Dots	12	HPV 39	23	HPV 56	34	HPV 70
2	Positive control	13	HPV 40	24	HPV 57	35	HPV 72
3	HPV 6	14	HPV 42	25	HPV 58	36	HPV 73
4	HPV 11	15	HPV 43	26	HPV 59	37	HPV 81 (CP8304)
5	HPV 16	16	HPV 44	27	HPV 61	38	HPV 82 (IS39)
6	HPV 18	17	HPV 45	28	HPV 62	39	HPV 82 (MM4)
7	HPV 26	18	HPV 51	29	HPV 66	40	HPV 83 (MM7)
8	HPV 31	19	HPV 52	30	HPV 67	41	HPV 84 (MM8)
9	HPV 33	20	HPV 53	31	HPV 68a	42	HPV 90
10	HPV 34	21	HPV 54	32	HPV 68b	43	HPV 91
11	HPV 35	22	HPV 55*	33	HPV 69		

*HPV 55 is classified by now as subtype of HPV44, but is still labelled HPV55 for consistency reasons.

	High Risk HPV		probably High Risk HPV
	Guide Dot, positive control		Low Risk HPV

In addition to the 41 HPV genotypes following controls are situated on the chip:

- Guide dots (Hybridization control): These dots are used by the *VisionArray* Analyzer Software for the positioning of the grid. Additionally, staining of the guide dots is proof for a successful hybridization, labeling, and staining reaction and is used for the calculation of the relative intensity of the signals.
- Positive control/PCR control: These controls are used for the evaluation of the PCR-reaction and the quality of the PCR-template.

If any of these controls fail to show complete staining, the results are not valid.

All capture sequences and the positive control are set up on the chip as duplicates and the guide dots as triplicates. The signals are visible on the chip as circular hybridization signals.

5 Interpretation of Results

5.1 General

With the aid of the *VisionArray* HPV Chip 1.0, it is possible to make a qualitative statement about the presence or absence of one or more of 41 HPV types in the investigated sample.

The intensity of the signals is influenced by the prevalence of the target sequences in the sample as well as miscellaneous factors of the detection system. The absolute numbers of the signal intensity cannot be used for quantification of the DNA-concentration.

5.2 Software-Based Evaluation

The automated evaluation of the results is performed by the *VisionArray* Analyzer Software. A comprehensive manual for a chip-analysis is enclosed to the Software.

5.3 Specificity of the Tests

5.3.1 Analytical Specificity

Analytical specificity of the VisionArray HPV Chip 1.0 was tested for each of the 41 HPV types separately. Therefore, sequence verified plasmids with a concentration of $\geq 100,000$ genome equivalents (GEQ) of the 40 other HPV types were tested in triplicates. For this purpose, all HPV types were amplified by PCR, hybridized on the chip and subsequently detected.

Specificity is calculated as followed:

$TN / (TN + FP)$ [*TN= True negative signals; FP= False positive signals*]

Specificity for all 41 HPV-types

HPV Type	Specificity [%]
6	100
11	100
16 (HR)	100
18 (HR)	100
26	100
31 (HR)	100
33 (HR)	100
34	100
35 (HR)	100
39 (HR)	100
40	99.2
42	100
43	100
44	100
45 (HR)	100
51 (HR)	100
52 (HR)	100
53	100
54	100
55	100
56 (HR)	100
57	100
58 (HR)	100
59 (HR)	100
61	100
62	100
66	100
67	100
68a	100
68b	97.6
69	100
70	100
72	100
73	100
81 (CP8304)	98.4
82 (IS39)	100
82 (MM4)	100
83 (MM7)	100
84 (MM8)	100
90	100
91	100

5.3.2 Cross hybridizations:

- In 33% of the cases HPV 70 hybridized to HPV 40. However, in lower concentrations no cross hybridization could be observed.
- In 66% of the cases HPV 62 hybridized to HPV 81. However, in lower concentrations no cross hybridization could be observed.
- In 100% of the cases HPV 68a hybridized to HPV 68b. In lower concentrations no cross hybridization could be observed. However, HPV 68b is a subtype and therefore in high due to the cross hybridisations not distinguishable from HPV 68a.

5.4 Sensitivity of the test

The sensitivity of the analysis system was tested for each HPV type separately. The sensitivity depends on the amount and efficiency of the PCR-cycles and the affinity of the catchers.

In order to estimate the sensitivity, cloned DNA-samples of the HPV types were tested in 3-6 different dilutions ranging from 50 GEQ to 500,000 GEQ. The lowest concentration that could be detected in at least 66 % of the cases was set as the limit of detection (LOD) of the analyte. Only for HPV 16 and 18 the WHO Human Papillomavirus Laboratory Manual (First edition, 2009) suggests a LOD of 50 GEQ.

Limit of Detection for all 41 HPV-types

HPV Type	Limit of Detection (GEQ)
6	50
11	50
16 (HR)	50
18 (HR)	50
26	500
31 (HR)	5,000
33 (HR)	50
34	50
35 (HR)	50
39 (HR)	50
40	50
42	500
43	500
44	500
45 (HR)	50
51 (HR)	50
52 (HR)	500
53	500
54	50
55	50,000
56 (HR)	50
57	500
58 (HR)	500
59 (HR)	5,000
61	500
62	500
66	500
67	50
68a	500,000
68b	500
69	500
70	50
72	5,000
73	5,000
81 (CP8304)	50
82 (IS39)	50
82 (MM4)	500
83 (MM7)	500
84 (MM8)	5,000
90	500
91	500

The determined sensitivity refers to the detection of a single target sequence. The detection of a multiple infection can lead to impairment of the sensitivity of some HPV types, due to competition during the PCR-reaction, especially in mixed probes with a strong difference in the concentration. In extremely high concentrated HPV-probes the intensity of the positive control might be impaired as well.

5.5 Cutoff

For the evaluation of the results the dot size is set to 50.

The threshold (cutoff) was set to 25 for the greyscale image of this dot size. A signal below this value is considered background by the *VisonArray* Analyzer software.

5.6 Limits of the Method

Interpretation of results must be made within the context of the patient's clinical history with respect to further clinical and pathologic data by a qualified pathologist.

Furthermore, following factors can influence the detection system:

- Deviation from the proposed detection protocol (e.g. temperature or volumes of the reagents)
- Degraded or low concentrated DNA material
- Inappropriate raw material
- Use of not calibrated or impaired equipment
- Low PCR efficiency due to PCR inhibition in DNA raw material (e.g. blood)
- Use of PCR additives that could influence the hybridization (e.g. DMSO, betaine, urea)

6 Literature

IARC Monographs on the evaluation of carcinogenic risks to humans, Vol. 100, 2012; ISBN 978 92 832 1319 2

WHO Human Papillomavirus Laboratory Manual, First edition, 2009.

Schmitt M, et al. (2008) *Journal of Clinical Microbiology* **46**:1050-1059.










Schmitt M, et al. (2013) *International Journal of Cancer* **132**:2395-2403.

7 Problems and Possible Causes

Any deviation from the operating instructions can lead to impairment of the detection reaction of the target sequence

Problem	Possible cause	Action
No signal	Wrong raw material	Check the raw materials
	Wrong combination of chip and sample	Check the sample/chip combination
	Wrong temperature	Check the hybridization temperature
	Expired reagents	Check the reagents
Only guide dots and no other signals	Problems with the PCR product (PCR not efficient enough or DNA template degraded)	Check PCR efficiency with a positive control; Check PCR chemicals and thermal cycler program; Check PCR product in agarose gel
Only guide dots and PCR control, but no other signals	No target sequence present	Use positive control
Too much background	Incubation time of Detection Solution or Blue Spot Solution too long; Temperature during incubation too high	Check incubation time and temperature of Detection Solution and Blue Spot Solution
Strong, leaking signals	Incubation time of Detection Solution or Blue Spot Solution too long	Stepwise adjustment of the incubation time of Detection Solution and Blue Spot Solution
Weak signals	Hybridization temperature incorrect	Check temperature
	Hybridization time too short	Extend hybridization time
	Incubation time of Detection Solution or Blue Spot Solution too short	Extend incubation time of Detection Solution and Blue Spot Solution
Cross-hybridization signals, false positive signals	Contamination of the PCR chemicals or PCR product	Replace the PCR chemicals in use
	Contamination during the preparation of the PCR or of the hybridization mix	Avoid transfer of sample during the preparation of the mix
	Hybridization temperature too low	Check hybridization temperature
	Several chips incubated too long in the same wash buffer	Swift execution of the washing steps
Single signal instead of duplicates	Mechanical elimination of the second signal, e.g. due to contact with the pipette tip	Avoid direct contact with the array field
	Irregular covering of the array field due to air bubbles	Apply solutions without air bubbles

8 Explanation of the Symbols

	Check information material
	Check Instruction for use
	Manufacturer
	Content sufficient for n tests
	LOT number
	Product code
	Use until
	Allowed storing temperature
	Store in the dark