

Heat Pretreatment Solution Citric

REF PT-0001-1000

Σ 14 (2x 500 ml)


For use in *in situ* hybridization procedures

CE

IVD

In vitro diagnostic medical device

according to EU directive 98/79/EC



1. Scope of Application

This product is designed for *in vitro* diagnostic use (according to EU directive 98/79/EC). Interpretation of results must be made within the context of the patient's clinical history with respect to further clinical and pathologic data of patient by a qualified pathologist!

The Heat Pretreatment Solution Citric (PT1) is designed to be used for heat pretreatment of formalin-fixed, paraffin-embedded tissue or cell samples in *in situ* hybridization (ISH) protocols of ZytoVision.

2. Safety Precautions and Disposal

- ✓ Read the operating instructions prior to use!
- ✓ Do not use the reagents after the expiry date has been reached!
- ✓ Avoid any cross-contamination and micro-bacterial contamination of the reagents!
- ✓ Avoid any direct contact with the reagents. Take appropriate protective measures (use disposable gloves, protective glasses, and lab garments)!
- ✓ If reagents come into contact with skin, rinse skin immediately with copious quantities of water!
- ✓ Never pipet solutions with your mouth!
- ✓ The disposal of reagents must be carried out in accordance with local regulations!
- ✓ A material safety data sheet is available on request for the professional user!

3. The Heat Pretreatment Solution Citric

The following components are included:

Code	Component	Quantity	Container
PT1	Heat Pretreatment Solution Citric	2x 500 ml	Screw-cap bottle (large)
	Instruction manual	1	

Component (**PT1**) is sufficient for 14 staining jars of 70 ml each.

4. Storage and Shelf Life

The component must be stored at 2...8°C. If these storage conditions are followed, the component will function, without loss of performance, at least until the expiry date printed on the label

5. Instructions

Heat Pretreatment Solution Citric (PT1) is a ready-to-use solution for heat pretreatment in ISH applications, following the ISH protocols of ZytoVision.

Our experts are available to answer your questions.

As of: June 22, 2010 (4.6)

Trademarks:

ZytoVision® is a trademark of ZytoVision GmbH.



ZytoVision GmbH · Fischkai 1
D - 27572 Bremerhaven · Germany
Phone: +49 (0)471/4832 - 300
Fax: +49 (0)471/4832 - 509
www.zytovision.com
info@zytovision.com



Your local distributor