




BIO SB TINTORETRIEVER PT MODULE

Epitope Retrieval System for use in Immunohistochemistry



*Ideal for automated workflow
environments*

*Intuitive touchscreen
interface*

Heat induced epitope retrieval and nucleic acid retrieval steps are most often performed using a conventional microwave, steamer, water bath or a pressure cooker. The Bio SB TintoRetriever PT is a precision-controlled heat source capable of maintaining a constant and reliable temperature while minimizing the potential for evaporation of the working solution under HIER procedure conditions. The TintoRetriever PT Module is recommended for use in both Immunohistochemistry (IHC), Immunocytochemistry (ICC) and In-Situ Hybridization (ISH) protocols.

Overview

- Intuitive and programmable touch interface.
- Slides can be quickly transferred from TintoRetriever PT Module to the TintoStainer.
- Generates epitope/nucleic acid retrieval report.
- Ideal for automated work flow environments.
- Rapid heat retrieval
- Intuitive touch screen interface
- Multiple temperature settings
- Use in IHC & ISH applications



| Product Description | Catalog # |
|-----------------------------------|-----------|
| TintoRetriever PT Module (110V) | BSB 7030 |
| TintoRetriever PT Module (220V) | BSB 7033 |
| Plastic Staining Dish | BSB 7009 |
| Slide Holder (24 Slides) | BSB 7010 |
| TintoRetriever Slide Thermometers | BSB 7005 |

TintoStainer Features



- Integrated bar-code scanning reduces user-error.
- Flexible protocol programming.
- Ideal for both research and clinical settings.
- Allows the use of antibody, reagent, and detection chemistry from any supplier.

| Product Description | Catalog # |
|--|-----------|
| TintoStainer 36 Slide (110v) with Report & Bar code Printers | BSB 7029 |
| TintoStainer 36 Slide (240v) with Report & Bar code Printers | BSB 7032 |
| TintoStainer 48 Slide (110v) with Report & Bar code Printers | BSB 7051 |
| TintoStainer 48 Slide (240v) with Report & Bar code Printers | BSB 7074 |
| TintoStainer 72 slide (110v) with Report & Bar code Printers | BSB 7052 |