

# BIO SB TINTORETRIEVER PRESSURE COOKER

Affordable, Precision Heat Epitope and Nucleic Acid Retrieval



**More Economical Than Other Suppliers!**

# TintoRetriever Pressure Cooker

The Bio SB TintoRetriever Pressure cooker is a precision-controlled heat source capable of maintaining a constant and reliable temperature while minimizing the potential for evaporation of the working solution. The TintoRetriever Pressure Cooker can be used for a variety of molecular pathology applications including IHC, FISH, and CISH.

Most formalin-fixed tissues require an antigen or nucleic acid retrieval step prior to any immunochemical staining. The retrieval process breaks some of the methylene bridges that crosslink proteins in formalin-fixed tissues and allow antibodies or RNA/DNA probes to bind to unmasked epitopes or nucleic acids. Heat epitope or nucleic acid retrieval are most often performed using a conventional microwave oven or water bath. The TintoRetriever Pressure Cooker is a solution for laboratories looking for a comprehensive, effective, and reliable retrieval system.

## Overview

- Cost effective.
- Fast, uniform and reproducible.
- Built in temperature gauge for validation.
- Supports CISH/FISH nucleic acid retrieval.
- Supports IHC, ICC and IF epitope retrieval.
- Handles up to 96 slides.
- Operates at 110 to 220 Volts.
- Accessories and replacement parts available!



**Above:** TintoRetriever Pressure Cooker Bundle, Includes 4 Plastic Staining Dishes, 4 Slide Holders, 1 Staining Dish Support, 200 ml ImmunoRetriever 20X Citrate, 200 ml ImmunoRetriever 20X EDTA.

| Product Description                    | Catalog # |
|--|-----------|
| TintoRetriever Digital Pressure Cooker | BSB 7087  |



**Left:** Staining Dish Support, holds 96 slides total!

| Product Description                          | Catalog # |
|--|-----------|
| Staining Dish Support                        | BSB 7086  |
| Plastic Staining Dish                        | BSB 7009  |
| Slide Holder, 24 places                      | BSB 7010  |
| TintoRetriever Slide Thermometers, Pack of 5 | BSB 7005  |